



The **Allen Yellofins Elite** stirrups enable safe and easy positioning while providing enhanced surgical access. Yellofin stirrups allow for easy adjustment in patients up to 500 lbs (227 kg). The boot on the Yellofin extends along its lateral side with a pad to protect the head of the fibula and peroneal nerve. The boot is self adjusting to minimize pressure on the calf. These stirrups feature lift assist technology.

## Features

- Extended lateral fin boot design.
- Self adjusting boot.
- Lift-assist technology.
- Squeeze grip handle provides easy intraoperative adjustment.
- Lithotomy and length indicators.



## Specifications

### General

**Patient Weight Capacity:** 500 lbs (227 kg)

**Lithotomy Range:** +84° to -33°

**Abduction Range:** +25° to -9°

### Procedures

GYN, Urological, Laproscopic, Robotic, and Colorectal procedures including:

**TAVH, LAVH**

**Colposcopy**

**Cryosurgery**

**LEEP Procedure**

**D&C**

**Radical Perineal Prostatectomy**

**General Surgery**

**AP Resections**

**Cysto Procedures**

**Radical Perineal Prostatectomy**

**Oncology Procedures**

**Robotic Surgical Procedures**

