



The **Allen Yellofins** stirrups enable surgical personnel to lift and position a patient's leg, then safely and securely hold the legs in position while surgical procedures are performed. The Yellofins have a lateral safety fin with patented lift-assist technology. The Design of the yellofin eliminates the possibility of peroneal nerve injury.

Features

- Squeeze grip handle provides easy intraoperative adjustment.
- Extended lateral fin boot design.
- Self adjusting boot.
- Lithotomy and length indicators.
- Lift-assist technology.



Specifications

General

Patient Weight Capacity: 350 lbs (159 kg)

Lithotomy Range: +84° to -33°

Abduction Range: +25° to -9°

Procedures

GYN, Urological, Laproscopic, Robotic, and Colorectal procedures including:

TAVH, LAVH

Colposcopy

Cryosurgery

LEEP Procedure

D&C

Radical Perineal Prostatectomy

General Surgery

AP Resections

Cysto Procedures

Radical Perineal Prostatectomy

Oncology Procedures

Robotic Surgical Procedures

