



The Axia Over Operating Table is a heavy-duty stainless steel, adjustable height table. It is manufactured with a fully welded construction with a reinforced worktop made with 304 series steel. The Over Operating table offers a manual height adjustment from 36.5 inches to 56 inches (depending on the model). The Table comes with 3-inch casters for easy movability.

Features

- Heavy duty 16 gauge #304 stainless steel construction
- 1.5" revealed edge for the strongest table
- 3" ball bearing casters (2 swivel, 2 locking)
- Adjustable height from 35" to 56" (depending on table Model)

Axia Over Operating Table Specifications

Models and Dimensions

- AX-HA203628 - 20"W x 36"L x 28" Leg Clearance
- AX-HA243628 - 24"W x 36"L x 28" Leg Clearance
- AX-HA244840 - 24"W x 48"L x 40" Leg Clearance
- AX-HA246052 - 24"W x 60"L x 52" Leg Clearance
- AX-HA463426 - 46"W x 34L x 26" Leg Clearance
- AX-HA3442 - 34"W x 42"L x 36.5" Leg Clearance

