



The Axia Closed Case Cart is a stainless steel welded cart that can securely and sterilely transport surgical instruments and supplies. The closed carts are built for strength and quality, offering a soundproof, quiet operation of the doors. This Axia cart comes in a variety of styles, including glass doors. The Glass in these doors is made of polycarbonate for high temperatures and high impact..

Features

- Fully welded construction with an inner skeletal system for strength and easy service with bolt-on features
- Soundproof tops and doors for quiet operation
- 6" maintenance-free, high temp, high-pressure stainless steel casters with sealed ball bearings (3 swivels & 1 fixed directional)
- Bolt-on exterior push handle for better control
- Bolt-on full wrap-around non-marring bumper assembly
- Bolt-on door hinges
- Bolt-on positive lock door latches
- Bolt-on 8.5" x 11" card holder
- Bolt-on tamper evidence security hasps



Axia Closed Case Cart Specifications

Optional Features

- Glass door(s)
- Etching for serial tag ID/tracking
- Etching for end user logo
- Adjustable solid shelf (non-roll out)
- Adjustable solid shelf (roll-out)
- Adjustable perforated shelf (non-roll out)
- Adjustable perforated shelf (roll-out)
- Adjustable wire shelf (non-roll out)
- Adjustable wire shelf (roll-out)
- Side extensions
- Flag system (one each of red, yellow and green)
- Clean/dirty indicator (slide or wheel)
- Total locking casters
- Tow hitch

Model Dimensions

- Model 1
 - Case - 25"D x 40.25"W x 32"H(63.5 cm D x 102.23 cm W x 81.28 cm H)
 - Overall - 26.875" D x 45.75" W x 39.625" H(68.262 cm D x 116.20 cm W x 100.65 cm H)
- Model 2
 - Case - 25"D x 40.25"W x 47"H(63.5 cm D x 102.23 cm W x 119.38 cm H)
 - Overall - 26.875" D x 45.75" W x 54.625" H(68.262 cm D x 116.20 cm W x 138.75 cm H)
- Model 3
 - Case - 25"D x 23.625"W x 32"H(63.5 cm D x 60 cm W x 81.28 cm H)
 - Overall - 26.875" D x 29.125" W x 39.625"H (68.262 cm D x 73.97 cm W x 100.64 cm H)
- Model 4
 - Case - 25"D x 23.625"W x 47"H (63.5 cm D x 60 cm W x 119.38 cm H)
 - Overall - 26.875" D x 29.125" W x 54.625"H (68.262 cm D x 73.97 cm W x 138.74 cm H)
- Model 5
 - Case - 25"D x36"W x 32"H (63.5 cm D x 91.44 cm W x 81.28 cm H)
 - Overall - 26.875" D x 41.5" W x 39.625" H (68.262 cm D x 105.41 cm W x 100.64 cm H)
- Model 6
 - Case - 25.875"D x 46.5"W x 32"H (65.72 cm D x 118.11 cm W x 81.28 cm H)
 - Overall - 27.75"D x 52"W x 39.625" H (70.48 cm D x 132.08 cm W x 100.64 cm H)

