



Cost-Effective Panel Further Expands DR Line-Up to Meet the Needs of Any Imaging Department. The Axia DR + Wireless Panel is a productivity solution providing exceptional image quality, high durability, and functionality. The 14x17 panel uses enhanced Gadolinium technology and has unique features like drop notification emails. The panels fit seamlessly into existing x-ray rooms of environments including hospitals, imaging centers and mobile imaging units, and further expands 's extensive line-up of imaging systems across the spectrum and vastly improve work-flow.



Specifications

Panel Specifications:

Weight 7.7lbs
Pixel Size 150um
Imaging Area 14"x17"
Dimensions 384x460x15.2mm
7.4 Million Pixels
Pixel Array 2304x2800
3.3lp/mm
A/D Conversion 14 Bits
Wireless connectivity
Cassette sized
Lithium Ion battery
Off-center Imaging
Carbon fiber construction

Software Specifications

Generator Integration
Patient Review
Reverse Black/White
Multi-Scale Contrast Enhancement and Dynamic Range
Optimized Image Processing
Parameters

