



Features

- Keyed temperature lock out
- Push button digital operation
- Large dual digital displays
- Displays actual temperature and set point at all times
- Separate digitally controlled heating chambers
- Rapid warm time (2-6 hours depending on load)
- Uniform heating
- Fully insulated
- Audible and Visual Over-Temp Alarm
- U.L./C.U.L. approved
- FDA Registered
- OSHPD Pre-Appro



Specifications

Overall Size

30"W x 26.5"D x 74.5"T
17.2 cubic feet of storage
3 adjustable shelves

Upper chamber I. D.

26"W x 23"D x 15.25"T
5.27 cubic feet of storage
Temperature range: 90-160 F.
Accuracy +/- 1 degree F. of set point

Lower Chamber I. D.

26"W x 23"D x 34.5"T
11.93 cubic feet of storage
Temperature range: 90-160 F.
Accuracy +/- 3 degree F. of set point

Electrical

120V, 12 amp., 50/60HZ
.75 kw/hr, 2557 btu/hr
UL/CUL listed
FDA registered