



The **Axia Functional Dual-Sided Portable Training Headwall** is designed for teaching labs and institutions. This portable headwall provides vacuum air, and simulated oxygen just by plugging the unit into a wall outlet. Providing an inexpensive way for simulation a patient room with medical gas services, while giving a movable solution.

## Features

- Provides vacuum, air, and electrical power in a portable solution.
- Locking casters, 2 fixed and 2 moveable.
- Easily moveable from one bed to the next.

## Options

### Add

- Integrated Non-Functional Nurse Call
- Integrated Non-Functional Code Blue
- Integrated White Dry Erase Board
- Oxygen Flowmeter
- Medical Air Flowmeter
- Vacuum Regulator [Continuous]
- Vacuum Regulator [Continuous/Intermittent]
- Vacuum canister and holder
- Additional Mounting Brackets as requested

### Delete

- Delete pump and/or compressor

