



The Axia Instrument Table is a stainless steel table with a fully welded construction and reinforced top. The Tabletop and under shelf are made with 300 series stainless steel. The Instrument table comes with 2-inch swivel casters and stainless steel legs that are 1.25-inches thick.

Features

- Fully welded construction
- Reinforced tops
- Stainless steel legs (1 1/4")
- Stainless steel top (type 300 series)
- Stainless steel under shelf (type 300 series)
- 2" swivel casters



Axia Closed Case Cart Specifications

Optional Features

- 3" swivel casters
- Locking casters
- Drawers under all sizes
- Guard rails on 3 sides (available on top or bottom shelf)
- Guard rails on all 4 sides (available on top or bottom shelf)
- H or U brace instead of under shelf
- Side extensions
- 1" lip on back
- Multiple under shelves

Model Dimensions

H-Brace

- AX-H1620 - 16"W x 20"L x 34"H (40.64 cm W x 50.8 cm L x 86.36 cm H)
- AX-H1630 - 16"H x 30"L x 34"H (40.64 cm W x 76.2 cm L x 86.36 cm H)
- AX-H1833 - 18"W x 33"L x 34"H (45.72 cm W x 83.82 cm L x 86.36 cm H)
- AX-H2048 - 20"W x 48"L x 34"H (50.8 cm W x 121.92 cm L x 86.36 cm H)
- AX-H2036 - 20"W x 36"L x 34"H (50.8 cm W x 91.44 cm L x 86.36 cm H)
- AX-H2020 - 20"W x 20"L x 34"H (50.8 cm W x 50.8 cm L x 86.36 cm H)
- AX-H2430 - 24"W x 30"L x 34"H (60.96 cm W x 76.2 cm L x 86.36 cm H)
- AX-H2436 - 24"W x 36"L x 34"H (60.96 cm W x 91.44 cm L x 86.36 cm H)
- AX-H2448 - 24"W x 48"L x 34"H (60.96 cm W x 121.92 cm L x 86.36 cm H)
- AX-H2460 - 24"W x 60"L x 34"H (60.96 cm W x 152.4 cm L x 86.36 cm H)
- AX-H2472 - 24"W x 72"L x 34"H (60.96 cm W x 182.88 cm L x 86.36 cm H)
- AX-H3036 - 30"W x 36"L x 34"H (76.2 cm W x 91.44 cm L x 86.36 cm H)
- AX-H3060 - 30"W x 60"L x 34"H (76.2 cm W x 152.4 cm L x 86.36 cm H)
- AX-H3072 - 30"W x 72"L x 34"H (76.2 cm W x 182.88 cm L x 86.36 cm H)
- AX-H3048 - 30"W x 48"L x 34"H (76.2 cm W x 121.92 cm L x 86.36 cm H)
- AX-S3420 - 2" Casters, 34"H x 20"D x 16"D (86.36 cm H x 50.8 cm D x 40.64 cm D)
- AX-S3448 - 2" casters, 34"H x 48"D x 24"D (86.36 cm H x 121.92 cm D x 60.96 cm D)



Model Dimensions

U-Brace Tables

- AH-U1620 - 16"W x 20"L x 34"H (40.64 cm W x 50.8 cm L x 86.36 cm H)
- AH-U1620 - 16"W x 30"L x 34"H (40.64 cm W x 76.2 cm L x 86.36 cm H)
- AH-U2436 - 24"W x 36"L x 34"H (60.96 cm W x 91.44 cm L x 86.36 cm H)
- AH-U2448 - 24"W x 48"L x 34"H (60.96 cm W x 121.92 cm L x 86.36 cm H)
- AH-2460 - 24"W x 60"L x 34"H (60.96 cm W x 152.44 cm L x 86.36 cm H)

MRI Tables

- AX-MRI1620 - 16"W x 20"L x 34"H (40.64 cm W x 50.8 cm L x 86.36 cm H) with shelf
- AX-MRI2020 - 20"W x 20"L x 34"H (50.8 cm W x 50.8 cm L x 86.36 cm H) with shelf
- AX-MRI1630 - 16"W x 30"L x 34"H (40.64 cm W x 76.2 cm L x 86.36 cm H) with shelf
- AX-MRI1833 - 18"W x 33"L x 34"H (45.72 cm W x 83.82 cm L x 86.36 cm H) with shelf
- AX-MRI2036 - 20"W x 36"L x 34"H (50.8 cm W x 91.44 cm L x 86.36 cm H) with shelf
- AX-MRI2048 - 20"W x 48"L x 34"H (50.8 cm W x 121.92 cm L x 86.36 cm H) with shelf
- AX-MRI2436 - 24"W x 36"L x 34"H (60.96 cm W x 91.44 cm L x 86.36 cm H) with shelf

