

New, Demo & Refurbished Medical Equipment

The **Axia Kite Series** surgical lights feature a touchscreen remote control, aluminum pedals, and adjustable petal movements. The Kite Series lights come in 6 different models. The Kite 4DD, Kite 4DS, Kite 3DM, Kite 3DD, Kite 3DS, Kite 3DM.

Features

- Aluminum petals.
- Adjustable petal movement.
- Touchscreen remote control.
- Intensity control.
- Autoclavable handle.
- Color temperature control.
- Stainless steel or painted body.

Phone: 1.800.GET.SOMA • www.SomaTechnology.com • Email: soma@somatechnology.com

Specifications

Kite 4DD



Intensity: 180,000 LUX

LED Average Life >50,000 Hrs.

Field Diameter: 100-150 mm

Number of LEDs: 88 x 2

Depth Illumination 100-150 mm

Intensity Controller: 10~100

Digital Color Temp.: 3,800-4,800 K

Power Supply: 100-240V 50/60Hz

CR Index: 93 Ra

Power Consumption: 44W x 2

Focus Adjustment: Manual

Kite 4DS



Intensity: 160,000 LUX

LED Average Life >50,000 Hrs.

Field Diameter: 100-150 mm

Number of LEDs: 88

Depth Illumination 100-150 mm

Intensity Controller: 10~100

Digital Color Temp.: 3,800-4,800 K

Power Supply: 100-240V 50/60Hz

CR Index: 93 Ra

Power Consumption: 44W

Focus Adjustment: Manual

Kite 4DM



Intensity: 160,000 LUX
LED Average Life >50,000 Hrs.
Field Diameter: 100-150 mm
Number of LEDs: 88
Depth Illumination 100-150 mm
Intensity Controller: 10~100
Digital Color Temp.: 3,800-4,800 K
Power Supply: 100-240V 50/60Hz
CR Index: 93 Ra
Power Consumption: 44W
Focus Adjustment: Manual



Specifications

Kite 3DD



Intensity: 160,000 LUX

LED Average Life >50,000 Hrs.

Field Diameter: 100-150 mm

Number of LEDs: 66 x 2

Depth Illumination 100-150 mm

Intensity Controller: 10~100

Digital Color Temp.: 3,800-4,800 K

Power Supply: 100-240V 50/60Hz

CR Index: 93 Ra

Power Consumption: 30W x 2

Focus Adjustment: Manual

Kite 3DS



Intensity: 120,000 LUX
LED Average Life >50,000 Hrs.
Field Diameter: 100-150 mm
Number of LEDs: 66
Depth Illumination 100-150 mm
Intensity Controller: 10~100
Digital Color Temp.: 3,800-4,800 K
Power Supply: 100-240V 50/60Hz
CR Index: 93 Ra
Power Consumption: 30W
Focus Adjustment: Manual

Kite 3DM



Intensity: 120,000 LUX

LED Average Life >50,000 Hrs.

Field Diameter: 100-150 mm

Number of LEDs: 66

Depth Illumination 100-150 mm

Intensity Controller: 10~100

Digital Color Temp.: 3,800-4,800 K

Power Supply: 100-240V 50/60Hz

CR Index: 93 Ra

Power Consumption: 30W

Focus Adjustment: Manual

