



## Features

- International LED Source, ultra-light, ultra-CRI, long life
- OPTical construction design, light spot, uniform
- All-in-one, ultra-light, high-end
- Specified product for surgeon doctors



## Specifications

### Dimensions

Height	250mm
Width	270mm
Depth	100mm
Weight	150g

### Power

Power	LED 3W
Power Supply Time	2.5h + 4.5h
Battery Charge Time	5h

### Light

Type	LED
Power	3W
Illumination	> 45,000lux
Color Temperature	5,500k
Color Rendering	85
Bulb Life	50,000h