



The **Axia Mayo Stands** are adjustable stainless steel stand with a removable tray top. This surgical mayo stand series can hold instruments and materials during the operating room or in-office procedures. The mayo stand height can easily be adjusted by either a knob or by foot.

Features

- Stainless steel and MRI compatible construction
- Knob or foot operation of adjustable height
- 2" swivel casters
- optional mobile base with 4 wheels

Models and Dimensions

Stainless Steel Mayo stands

- AX-MS24MB - 19.125" x 12.625" tray, 24" diameter, Mobile Base, Manual knob operated
- AX-MS16MB - 19.125" x 12.625" tray, 16" diameter, Mobile Base, Manual knob operated
- AX-MS730 - 19.125" x 12.625" tray, Manual knob operated
- AX-MS760 - 19.125" x 12.625" tray, foot operated

MRI Copatable Mayo Stands

- AX-MS-3-MRI - 19.125" x 12.625" tray, Manual knob operated
- AX-MS-3-MRI-MB - 19.125" x 12.625" tray, Manual knob operated, Mobile base

