

Medical / Crash Carts

# AXIA MEDICAL CARTS

## ABOUT AXIA CARTS

Axia Medical Crash Carts are a great solution for high-demand medical facilities. They offer additional space, optimal organization, and functionality. They possess from 5 to 9 drawers, and many other features depending on the model you pick.

## ABOUT AMC60R/50R CRASH CART

The AMC60R and AMC50R crash carts feature a durable construction of 6 drawers or 5 drawers allowing for an abundance of storage.

## ABOUT AMC60B ANESTHESIA CART

The AMC60R crash cart is an anesthesia cart which features a durable construction of 6 drawers allowing for an abundance of storage.

## ABOUT AMC90C PEDIATRIC CART

The AMC60R crash cart is a pediatric cart which features a durable construction of 9 drawers allowing for an abundance of storage.

## FEATURES

- Extension table
- 4 hook IV Pole
- Drawer dividers
- Seal lock or padlock option
- Noise-resistant castors with brakes
- Comfortable push handle
- Easy glide drawers
- Complete aluminum construction
- CPR board
- Oxygen tank holder
- External power strip





Axia Medical Carts

# SPECIFICATIONS

## AMC60R SPECIFICATIONS

### DIMENSIONS

Height: 44 in (111.8 cm)

Width: 26.4 in (67.1 cm)

Depth: 19.2 in (48.7 cm)

### DRAWERS

Drawers: 6

### DRAWER HEIGHTS

3 Drawers: 4.3 in (10.9 cm)

2 Drawers: 5.3 in (13.5 cm)

1 Drawer: 7.8 in (19.8 cm)

## AMC60B SPECIFICATIONS

### DIMENSIONS

Height: 44 in (111.8 cm)

Width: 26.4 in (67.1 cm)

Depth: 19.2 in (48.7 cm)

### DRAWERS

Drawers: 6

### DRAWER HEIGHTS

3 Drawers: 4.3 in (10.9 cm)

2 Drawers: 5.3 in (13.5 cm)

1 Drawer: 7.8 in (19.8 cm)

## AMC50R SPECIFICATIONS

### DIMENSIONS

Height: 39 in (99.1 cm)

Width: 26.4 in (67.1 cm)

Depth: 19.2 in (48.7 cm)

### DRAWERS

Drawers: 5

### DRAWER HEIGHTS

3 Drawers: 4.3 in (10.9 cm)

1 Drawers: 5.3 in (13.5 cm)

1 Drawer: 7.8 in (19.8 cm)

## AMC90C SPECIFICATIONS

### DIMENSIONS

Height: 46 in (111.8 cm)

Width: 26.4 in (67.1 cm)

Depth: 19.2 in (48.7 cm)

### DRAWERS

Drawers: 9

### DRAWER HEIGHTS

8 Drawers: 3.75 in (9.5 cm)

1 Drawer: 7.8 in (19.8 cm)