



The Axia Nesting Table is a stainless steel surgical table. The Nesting table is built with a fully welded construction. The table top is 304 series stainless steel. The table legs and U brace stabilize the table so that it is easy to move. It comes with 2-inch, non-locking, swivel casters.

Features

- Fully welded construction
- Stainless steel legs (1 1/4")
- Stainless steel tops (type 304)
- Stainless steel U brace (1 1/4")
- 2" swivel casters

Axia Nesting Table Specifications

Model Dimensions

- AX-N95: 14"D x 28"W x 32"H
(35.56 cm D x 71.12 cm W x 81.28 cm H)
- AX-N94: 16"D x 32"W x 34"H
(40.64 cm D x 81.28 cm W x 86.36 cm H)
- AX-N93: 18"D x 36" W x 36"H
(45.72 cm D x 91.44 cm W x 91.44 cm H)
- AX-N92: 20"D x 40"W x 38"h
(50.8 cm D x 101.6 cm W x 96.52 cm H)
- AX-N91: 22"D x 44"W x 40"H
(55.88 cm D x 111.76 cm W x 101.6 cm H)
- AX-N90: 24"D x 48"W x 42"H
(60.96 cm D x 121.92 cm W x 106.68 cm H)
- AX-N69: Nesting Table,
(Set Of 6 Nested Tables)

