



Axia P3 Podiatry Exam Chair - Power up/ down, front/ back tilt, reclining up/ down. Pneumatic assist single leg support. 180 degree chair rotation. Foot control. Available in grey, blue, and more to come.

## Specifications

<b>Actuators</b>	3
<b>Length</b>	70" - 78"
<b>Chair width</b>	32"
<b>Sitting width</b>	24"
<b>Height</b>	19-42"
<b>Chair Weight</b>	63lbs
<b>Max. Capacity</b>	440 lbs
<b>Max. Working Height</b>	52"
<b>Watt</b>	96
<b>Fold Away Arm Support</b>	185°
<b>Memory button</b>	No
<b>Rotatable</b>	Yes

