



The **Axia SB1** imaging table is cantilevered (diving board) style design. The “H” style base allows maximum C-Arm clearance along the entire cantilevered length. It comes standard with 500 pound load rating, 2” patient pad, 24” wide top, locking swivel castors, narrow base for easy patient transfer, low boarding height.

Features

- Cantilevered (diving board) style design.
- “H” style base.
- 500 pound load rating.
- 2” patient pad.
- 24” wide top.
- Low boarding height.



Specifications

Dimensions

Height: 28 - 38" (71.12 - 96.52 cm)

Width: 26" (66.04 cm)

Length: 86" (218.44 cm)

Weight: 520 lbs (236 kg)

Tabletop

Height: 1.62" (4.11 cm)

Width: 24" (60.96 cm)

Length: 84" (213.36 cm)

Mattress: 2" (5.08 cm)

Tabletop and Mattress attenuation: 1.6mm Al.(+/- .2mm tolerance)

Material: Carbon fiber skin/foam

Patient Straps: 2

Radiolucent Area

Width: 24" (60.96 cm)

Length: 56" (142.24 cm)

Steel Base

Shape: H-Shaped

Length: 51" (129.54 cm)

Width: 26" (66.04 cm)

Finish: Epoxy

Castors: 5" full swivel and wheel lock

Movements

Height: 28 - 38" (71.12 - 96.52 cm)

Hand Control: Activates Motion

