





The **Axia SL2** imaging table is designed for high patient throughput for surgery centers, pain management, and other multidisciplinary medical facilities. It features a diving board style for easy C-arm access and a large patient capacity of 500 pounds. The Axia SL2 is a 2 movement table. This table is available with 2 different style table tops; standard and cervical. The cervical top is 22 inches wide with a fully adjustable headrest.

Features

- Patient Capacity: 500lb (227kg)
- Tabletop Material: Carbon fiber skin/foam
- core with integral facial cut-out
- Radiolucent Area: 56" x 24"
- Tabletop: 84"L x 24"W
- Table Dimensions: 86"L x 26"W
- 5" Castors
- Mattress: 2" thick (5 cm)
- Available with a standard or cervical top.



Standard Top



Cervical Top

Specifications

Dimensions

Table Length: 86 in (218.4 cm)

Table Width: 26 in (66 cm)

Table Height: 28 - 38 in (71.1-96.5 cm)

Weight Capacity: 500 lbs (226.8 kg)

Tabletop / Imaging Area / General

Tabletop (LWH): 84 x 24 x 1.62 in (213.4 x 61 x 4.1 cm)

Radiolucent Area (LW): 56 x 24 in (142.2 x 61 cm)

Material: Carbon Fiber - skin/foam

Casters: 5 in (12.7 cm)

Base (LW): I-Shaped Steel - 82 x 26 (208.3 x 66 cm)

Hand Control: Activates all Motions

Patient Straps: 2

Mattress: 2 in (5.1 cm)

Facial Cutout: Yes

Power

AC Voltage: 120V~

Frequency: 60 Hz

AC Current: 2.5 A

Electric Shock Protection: Class I

Applied Part Rating: Type B

Mode of Operation: 10% Duty Cycle; 1 min on, 9 minutes off

Movements

Height: 28 - 38 in (71.1-96.5 cm)

Lateral Tilt: $\pm 17^\circ$

