





The **Axia STL4** is a mobile imaging table which is intended to support a 500 pound patient. The Axia STL4 is sold new and is intended to be used with mobile x-ray systems for minimally invasive procedures. The STL4 is made in the USA and has 4 table movements. This table is available with 2 different style table tops; standard and cervical. The cervical top is 22 inches wide with a fully adjustable headrest.

## Features

- Easy c-arm access with diving board style table top.
- Large patient capacity up to 500 pounds.
- 2-inch thick mattress pad.
- Facial cutout for prone procedures.
- Easy patient transfer with low height.
- Available with a standard or cervical top.



**Standard Top**



**Cervical Top**

# Specifications

**Tabletop Dimensions**      **Length:** 84" (214.4 cm)  
    **Width:** 24" (61cm)  
    **Thickness:** 1.62" (4.1 cm)

**Radiolucent Dimensions**      **Length:** 84" (142.2 cm)  
    **Width:** 24" (61cm)

**Table Dimensions**      **Length:** 86" (218.4 cm)  
    **Width:** 26" (66 cm)

**Weight**      **Table Weight:** 590lbs (268kg)  
    **Patient Capacity:** 500lbs (227kg)

**Table Movements**      **Height:** 30" to 40" (76.2 to 101.6 cm)  
    **Lateral Tilt:** ( + / - ) 17°  
    **Trendelenburg:** 13°  
    **Longitudinal Travel:** 10" (+/- 5) (25.4 cm (+/- 12.7 cm))

**General**      **Mattress Thickness:** 2" (5 cm)  
    **Tabletop / Mattress Attenuation:** 1.6mm Al.(+/- .2mm tolerance)  
    **Tabletop Material:** Carbon fiber skin/foam core with integral facial cut-out  
    **Casters:** 5" casters with full swivel and wheel lock  
    **Steel Base:** I-shaped with epoxy finish  
    **Hand Control:** Activates all motions

**Power**      120V/60HZ or 230V/50HZ  
    Protect product from moisture  
    Allowable Storage & Operating Humidity: 0-92%  
    Allowable Storage Temperature: 32°F to 104°F (0°C to 40°C)  
    Allowable Operating Temperature: up to 104°F (40°C)

