


Features

-



Specifications

Dimensions

Height: 33.5 to 43.5" (85.1 to 110.5 cm)

Width: 26" (66 cm)

Length: 86" (218.4 cm)

Weight: 695 lbs (315 kg)

Tabletop

Material: Carbon fiber skin/foam core with integral facial cut-out

Mattress Thickness: 2" (5.1 cm)

Length: 84" (213.4 cm)

Width: 24" (61 cm)

Patient Straps: 2

Radiolucent Area

Length: 56" (142.2 cm)

Width: 24" (61 cm)

Base

Shape: I Shaped

Length: 82" (208.3 cm)

Width: 26" (66 cm)

Castors: 5" with full swivel and wheel lock

Finish: Epoxy Finish

Movements

Height: 33.5" to 43.5" (85.1 to 110.5 cm)

Lateral Tilt: $\pm 17^\circ$

Trendelenburg: $\pm 13^\circ$

Longitudinal travel: $\pm 5''$ (± 12.7 cm)

Lateral travel: Total: 10" ($\pm 5''$); Total: 25.4 cm (± 12.7 cm)

