



Features

- Keyed temperature lock out
- Push button digital operation
- Large dual digital displays
- Displays actual temperature and set point at all times
- Separate digitally controlled heating chambers
- Rapid warm time (2-6 hours depending on load)
- Uniform heating
- Fully insulated
- Audible and Visual Over-Temp Alarm
- U.L./C.U.L. approved
- FDA Registered
- OSHPD Pre-Appro



Specifications

Overall Size

Size: 24"W x 20.5"D x 24.5"T

Chamber I. D.

20"W x 17"D x 15.25"T
3.0 cubic feet of storage
Temperature range: 90-160 F.
Accuracy +/- 1 degree F. of set point

Electrical

120V, 7.5 amp., 50/60HZ
.45 kw/hr, 1537 btu/hr
UL/CUL listed
FDA registered

