



The Axia Sealed Case Cart is a stainless steel surgical cart that can transport surgical instruments and supplies around your facility. The sealed cart is manufactured with a fully welded construction that offers a soundproof top and doors for quiet operation. The cart has two external push handles for better control when transporting.

Features

- Fully welded construction
- Soundproof tops and doors for quiet operation
- Exterior push handle for better control
- Door gaskets for increased dust control
- Locking rods to keep doors secure
- Full wrap around non-marring heavy-duty bumper
- Magnetic door catches to keep doors open while in use
- 8.5" x 11" card holder
- Tamper evidence security tabs
- 6" maintenance-free, high-temperature, high-pressure stainless steel casters with sealed ball bearings (2 rigid & 2 swivel with brake)



Axia Sealed Case Cart Specifications

Optional Features

- Etching for serial tag ID/tracking
- Etching for end user logo
- Adjustable solid shelf (non-roll out)
- Adjustable solid shelf (roll-out)
- Adjustable perforated shelf (non-roll out)
- Adjustable perforated shelf (roll-out)
- Adjustable wire shelf (non-roll out)
- Adjustable wire shelf (roll-out)
- Side extensions
- Flag system (one each of red, yellow, and green)
- Clean/dirty indicator (slide or wheel)
- Total locking casters
- Tow hitch

Model Dimensions

- Model 1
 - Case - 25.5"D x 41.25"W x 32"H (64.77 cm D x 104.775 cm W x 81.28 cm H)
 - Overall - 27.375"D x 46.75"W x 39.5"H (69.532 cm D x 118.745 cm W x 100.33 cm H)
- Model 2
 - Case - 25.5"D x 41.25"W x 47"H (64.77 cm D x 104.775 cm W x 119.38 cm H)
 - Overall - 27.375"D x 46.75"W x 54.5"H (69.532 cm D x 118.745 cm W x 138.43 cm H)
- Model 3
 - Case - 25.5"D x 24.625"W x 32"H (64.77 cm D x 62.547 cm W x 81.28 cm H)
 - Overall - 27.375"D x 30.125"W x 39.5"H (69.532 cm D x 76.517 cm W x 100.33 cm H)
- Model 4
 - Case - 25.5"D x 24.625"W x 47"H (64.77cm D x 62.547 cm W x 119.38 cm H)
 - Overall - 27.375"D x 30.125"W x 54.5"H (69.532 cm D x 76.517 cm W x 138.43 cm H)
- Model 5
 - Case - 25.5"D x 37"W x 32"H (64.77 cm D x 93.88 cm W x 81.28 cm H)
 - Overall - 27.375"D x 42.5"W x 39.5"H (69.53 cm D x 107.95 cm W x 100.33 cm H)
- Model 6
 - Case - 26.375"D x 47.5"W x 32"H (66.99 cm D x 120.65 cm W x 81.28 cm H)
 - Overall - 28.25"D x 53"W x 39.5"H (71.755 cm D x 134.62 cm W x 100.33 cm H)

