



The Axia Space Saver Table is a stainless steel welded table that can hold surgical instruments and supplies. The table offers options for a fold-down top shelf and a shelf underneath the table. It comes with lock and swivel casters for easy movability.

Features

- Stainless steel tilted upper shelf to easily view instruments
- Non-adjustable height and adjustable height models
- 4" Casters, 2 locking, 2 swivel
- 30" depth table top

Space Saving Table Specifications

Non-Adjustable Height Tables

- AX-TT3048 - 30"D x 48"W x 34"H
(76.2 cm D x 121.92 cm W x 86.36 cm H)
- AX-TT3060 - 30"D x 60"W x 34"H
(76.2 cm D x 152.4 cm W x 86.36 cm H)
- AX-TT3072 - 30"D x 72"W x 34"H
(76.2 cm D x 182.88 cm W x 86.36 cm H)
- AX-TTSS603030"D x 60"D x 36"H
(76.2 cm D x 152.4 cm W x 91.44 cm H)
- AX-TTSS723030"D x 72"W x 36"H
(76.2 cm D x 182.88 cm W x 91.44 cm H)

Adjustable Height Tables

- Ax-TTSS6030-EH - 30"D x 60"W x 36"H
(76.2 cm D x 152.4 cm W x 91.44 cm H)
Electronic adjustable Height

