





The **Axia TL3** is a mobile imaging table which is intended to support a 500 pound patient. The Axia TL3 is intended to be used with mobile x-ray for minimally invasive procedures. This table is available with 2 different style table tops; standard and cervical. The cervical top is 22 inches wide with a fully adjustable headrest.

Features

- The Axia TL3 is a **NEW** mobile imaging table which is intended to support a 500 pound patient.
- The Axia TL3 is intended to be used with mobile x-ray for minimally invasive procedures.
- 3 Move Imaging Table:
 - Motor Driven 28" to 38" Height Adjustment
 - Motor Driven Lateral Tilt: $\pm 17^\circ$
 - Motor Driven Trendelenburg: $\pm 13^\circ$
- Made in the USA
- Available with a standard or cervical top



Standard Top



Cervical Top

Specifications

Dimensions

Table Top: 84"L x 24"W x 1.62"T

Table Dimensions: 86"L x 26"W

General Specifications

Hi/Low

Trendelenburg/Reverse Trendelenburg

Lateral Tilt

24" x 84" Carbon Fiber Diving Board

500 lb. Patient Weight Capacity

5" Locking Casters

Hand Held Remote Control

120VAC, 60Hz., 2.5 Amp Power

