



Axia's **UroMax 4** C-arm surgical tables offer an excellent solution for real-time urological X-ray procedures. The cantilevered, carbon fiber, radiolucent top areas allow full fluoroscopic visualization and unobstructed C-arm positioning. The UroMax 4 has highly versatile movements, options, and accessories. The carbon fiber leg extension accessory can quickly turn your Urology table into a procedure table for pain management and other C-arm guided procedures.

Features

- Carbon fiber with integral facial cutout section
- Leg extension easy to remove and very light weight
- Diving board style for easy C-arm access
- Large patient capacity (500-lb pound load)
- A 2" patient high-quality comfort pad
- Light weight – easy to physically move the table
- Our 24" wide top eliminates the needs for arm boards
- Easy patient transfer to table
- Low height for patient boarding

*Stirrups are an optional accessory



Specifications

Dimensions

Table Length: 55 in (139 cm)
Table Width: 24 in (60.9 cm)
Table Height: 30 – 40 in (76.2 – 101.6 cm)
Weight Capacity: 500 lbs (226.7 kg)
Table Weight: 740 lbs (335.7 kg)
H-Shaped Steel Base: 51 x 26 in (129.5 – 66 cm)
Mattress: 2 in thick (5 cm)
Radiolucent Area: 27 x 24 in (68.5 x 60.9 cm)
Table Top: 55 x 24 x 1.62 in (139 x 60 x 4 cm)

Tilts

Lateral Tilt: $\pm 17^\circ$
Trendelenburg: $\pm 13^\circ$
Longitudinal Travel: 10 in (± 5 in)

Power

AC Voltage: 120V~
Frequency: 60 Hz
AC Current: 6.5A
Electric Shock Protection: Class I
Applied Part Rating: Type B
Mode of Operation: 10% Duty Cycle; 1 min on, 9 min off

Environmental

Storage Humidity: 0- 92%
Operating Humidity: 0- 92%
Storage Temperature: 32°F – 104°F (0°C – 40°C)
Operating Temperature: Up to 104°F (40°C)

