



The **Axia UroMax 5** C-arm surgical tables offer an excellent solution for real-time urological X-ray procedures. The cantilevered, carbon fiber, radiolucent top areas allow full fluoroscopic visualization and unobstructed C-arm positioning. The Axia UroMax 5 has highly versatile movements, options, and accessories. The carbon fiber leg extension accessory can quickly turn your Urology table into a procedure table for pain management and other C-arm guided procedures.

## Features

- Carbon fiber with integral facial cutout section
- Leg extension easy to remove and very lightweight
- Diving board style for easy C-arm access
- Large patient capacity (500-lb pound load)
- A 2" patient high-quality comfort pad
- Lightweight – easy to physically move the table
- Our 24" wide top eliminates the needs for arm boards
- Easy patient transfer to table
- Low height for patient boarding



\*Stirrups are an optional accessory

# Specifications

## Dimensions

**Table Length:** 55 in (139 cm)  
**Table Width:** 24 in (60.9 cm)  
**Table Height:** 32 – 40 in (81.2 – 101.6 cm)  
**Weight Capacity:** 500 lbs (226.7 kg)  
**Table Weight:** 740 lbs (335.7 kg)  
**H-Shaped Steel Base:** 51 x 26 in (129.5 – 66 cm)  
**Mattress:** 2 in thick (5 cm)  
**Radiolucent Area:** 27 x 24 in (68.5 x 60.9 cm)  
**Table Top:** 55 x 24 x 1.62 in (139 x 60 x 4 cm)

## Tilts

**Lateral Tilt:**  $\pm 17^\circ$   
**Trendelenburg:**  $\pm 13^\circ$   
**Longitudinal Travel:** 10 in ( $\pm 5$  in)  
**Latitudinal Travel:** 10 in ( $\pm 5$  in)

## Power

**AC Voltage:** 120V~  
**Frequency:** 60 Hz  
**AC Current:** 6.5A  
**Electric Shock Protection:** Class I  
**Applied Part Rating:** Type B  
**Mode of Operation:** 10% Duty Cycle; 1 min on, 9 min off

## Environmental

**Storage Humidity:** 0- 92%  
**Operating Humidity:** 0- 92%  
**Storage Temperature:** 32°F – 104°F (0°C – 40°C)  
**Operating Temperature:** Up to 104°F (40°C)

