




The **Axia V-MAX** is a diving board style four-way float imaging table for easy C-arm access. The Axia V-Max has a patient weight capacity of 500 pounds and a 2-inch high-quality pad. Different options include adding additional side rail systems, arm boards, or a foot switch. The table top features an instant response electric brake which releases for the four-way float function to operate.

Features

- Facial cut-out access.
- 24-inch wide 84-inch long carbon fiber table top.
- Joystick hand control.
- Rotating lockable casters.
- Longitudinal Travel: 36-inches
- Lateral Travel: 14-inches



Specifications

Dimensions

Height: 28-38 in (71.1-96.5 cm)

Width: 24 in (61 cm) (*Optional 27 in wide top (68.6 cm)*)

Length: 84 in (213.4 cm)

Weight Capacity: 500 lbs. (226.8 kg)

Movements

Height: 28-38 in (71.1-96.5 cm)

Longitudinal Travel: 36 in (91.4 cm)

Lateral Travel: 14 in (35.6 cm)

