



The Axia Vertical Handle Case Cart is a stainless steel cart that can transport surgical instruments and supplies. The Vertical Handles on each corner of the cart make it easier to steer and control in tight spaces. The cart offers a soundproof top and quiet operation doors. This surgical case cart also has a bumper that wraps around the entire cart.

## Features

- Ergonomically friendly vertical handles for precise control
- Fully welded construction
- 8.5" x 11" card holder
- Sound-proofed tops and doors for quiet operation
- Full wrap around non-marring bumper
- Tamper evidence security tab
- Positive lock door latches
- 6" maintenance-free, high temp, high-pressure stainless steel casters with sealed ball bearings (3 swivel & 1 directional lock)



# Axia Vertical Handle Case Cart Specifications

## Optional Features

- 6" Stainless steel total locking caster (1 swivel, 2 total lockings, 1 directional)
- Washable cart, high temp
- Cool touch handles (set of handles)
- Solid, adjustable roll-out shelf
- Perforated, adjustable non-roll out shelf
- Fixed top shelf
- Silicone door gasket

## Model Dimensions

- **AXVH-1072PS** - 26"D X 48.75"W X 51"H (66.04 cm D x 123.825 cm W x 129.54 cm H) / 24.75"D X 40"W x 38.75"H (62.865 cm D x 101.5 cm W x 98.425 cm H), 2PS, 2 Door
- **AXVH-952PS** - 27.75"D X 48.25"W X 53.5"H (70.485 cm D x 122.55 cm W x 135.89 cm H) / 26.5"D X 39.625"W X 40.75"H (67.31 cm D x 100.64 cm W x 103.50 cm H), 2PS, 2 Door

