



The Axia Work Table is a stainless steel surgical work table. The Work table is manufactured with a fully welded construction with a reinforced worktop made with 300 series steel. The 2-inch square steel legs and under-shelf add additional support to the table. The work table comes with 5-inch casters for easy movability.

## Features

- Fully welded construction
- Reinforced stainless steel top (type 300 series)
- Stainless steel under shelf (type 300 series)
- 2" square legs
- Added supports
- 5" casters (2 swivel & 2 locking)



# Axia Work Table Specifications

## Optional Features

- Ergonomically friendly adjustable heights (manual or electronic)
  - Manually adjustable height range: 36"-55"
  - Electronically adjustable height range: 36"-43"
- H or U brace instead of under shelf
- Full base unit underneath
- Drawers under all sizes
- Stainless steel guard rails
- Multiple under shelves

## Model Dimensions

- Model 1: 24"D x 36"W x 34"H (60.96 cm D x 91.44 cm W x 86.36 cm H)
- Model 2: 24" D x 48"W x 34"H (60.96 cm H x 121.92 cm W x 86.36 cm H)
- Model 3: 24:D x 60"W x 34"H (60.96 mc D x 152.4 cm W x 86.36 cm H)
- Model 4: 30"D x 72"W x 34"H (76.2 cm D x 182.88 cm W x 86.36 cm H)
- Model 5: 36"D x 72"W x 34"H (91.44 cm D x 182.88 cm W x 86.36 cm H)
- Modle 6: 0"D x 60"W x 34"H (76.2 cm D x 152.4 cm W x 86.36 H)

