



Designed for acute care adult and pediatric facilities, the compact 9" x 6" Infusomat® Space large volume pump weighs just 3 pounds. The lightweight large volume pump allows for flexibility and easy transport within a number of care settings. Along with industry-leading safety features, the lightweight mobile design allows you to quickly and easily stack up to 24 pumps in a single bedside tower. Equipped with an extensive drug library, easy navigation and intuitive interface, this small infusion pump packs state-of-the-art versatility that promotes efficiencies across numerous care settings.

Features

- Easy-to-learn, easy-to-use intuitive user interface
- 1,500-entry drug library available
- Alarms are audible, clearly visible and differentiated using different tones
- Available with or without outbound wireless capabilities
- Strong lightweight pole clamp holds up to three pumps mounted horizontally or vertically



Specifications

Dimensions	214x68x124 mm (WxHxD)
Weight	approx. 1.4 kg
Battery Type	Rechargeable NiMH battery pack, quick and easy exchange without opening device
Battery operating time	Approx. 4 h at 100 ml/h. When prompted, display shows battery status in h and min

