



Bair Hugger therapy temperature management units have become the standard of excellence in forced-air warming. The Bair Hugger Model 500/OR provides safe, quiet and effective warming to patients around the world. The temperature management units are the **ONLY** forced-air warming systems that can be used with an integrated, cost-effective method of fluid warming.

Features

- Small size frees up more usable workspace
- Light weight makes transport and set-up a breeze
- Freestanding or easily attached to an IV pole, bedrail or optional rolling cart
- Built-in hour meter makes it easy to monitor usage for preventative maintenance
- Unique snap-fit hose swivels at three points for easy blanket attachment and positioning



Specifications

Dimensions

Height: 24 in (61 cm)
Width: 14 in (36 cm)
Depth: 16 in (41 cm)
Weight: 43 lbs (19.5 kg)

Relative Noise Level

52 Decibels

Hose

Detachable, Flexible

Filtration System

0.2 M Levels

**Recommended Filter
Change**

Every 6 months or 500 hours

Device Ratings

110-120 VAC, 60 Hz, 9.5 Amperes, or
220-240 VAC, 50 Hz, 4.5 Amperes, or
100 VAC, 50/60 Hz, 10 Amperes

