

Baxa Microfuse Dual Rate Infuser Infusion/Syringe Pump



The **Baxa Microfuse Dual Rate Infuser** is intended for small volume infusion (SVI) of intravenous fluids. It is designed to deliver fluids that are diluted in 1 to 100 mL volumes and given intermittently over 1 to 120 minutes. Syringe infusion provides an alternative to gravity and minibag IV systems. The Microfuse Infuser deliver reliable, fixed-rate infusion using less costly materials, reducing waste and setup and administration labor.

Features

- Durable solid body construction with secure motor attachment.
- Easy-to-use simple push-button controls.
- Portable, lightweight, and slender profile.
- Reliable and low maintenance, with built-in safety features.
- Operates on 2 standard alkaline C-cell batteries.
- Average life of 400 hours.



