



The **BioMed Crossvent 4 Plus** offers ventilation support for patients ranging from infant, pediatric and adult. The ventilator comes with three different configuration settings; A standard setting, a setting with an air/oxygen blender attached, and one with an Air Entrainment. The Setting that offers an Air Entrainment setting can deliver 100 % more oxygen than most normal settings that can only deliver up to 50% without an air supply source. Using the color touch screen display, the operator can easily select a ventilation mode along with setting the various alarm modes. The Alarm modes include; tidal volume, PIP, PEEP, Mean pressure, low supply pressure, rate, and oxygen concentration. The Crossvent 4+ can be used in setting like patient transport, emergency room, intensive care, and many more respiratory needs.

Features

- Portable ventilator
- Infant-Adult
- SIMV
- CPAP
- PS
- Universal
- Color Display
- Comes either standard blender or with air entrainment

