

Buffalo Filter PlumeSafe Whisper Turbo Surgical Smoke Evacuators



The **Buffalo Filter PlumeSafe Whisper Turbo** is designed to swiftly and efficiently remove surgical smoke and noxious odors during laser and electrosurgical procedures. Employing advanced filtration technology, this system effectively captures and contains harmful particles, preventing their dispersion into the surgical environment. Its versatility ensures compatibility with a wide range of surgical applications, while its compact, portable design facilitates easy integration into any healthcare setting. With an intuitive control interface, it offers user-friendly operation and low noise output for a serene surgical atmosphere.

Features

- 4-Stage filtration in one casing.
- User-friendly LED control display.
- Turbo Mode - Significantly increases the air flow to quickly remove the smoke.
- Adjustable Flow Control - Suction setting can easily be adjusted in 10% increments.
- Filter housing and handles remain temperate with the outer surface to prevent burns.
- Quiet operation.



Buffalo Filter PlumeSafe Whisper Turbo Surgical Smoke Evacuators

Specifications

Dimensions

Height: 13.5 in (34 cm)

Width: 8 in (20 cm)

Depth: 17 in (43 cm)

Weight: 19.5 lbs (8.9 kg)

Display

Type: LED

Displays: Filter Status, flow rate, footswitch control

Filter Life Time

Factory set 6 hours (When using a 12 hour ViroSafe EY filter, the filter life timer will need to be reset for the second 6 hours of filter life.)

Variable Control

Yes

Static Suction

85 Inches H2O (211 mbar)

Noise Level

Maximum: 68 dB

