



The **SchureMed Candy Cane Stirrups** help position patients in a cost-effective way for a wide-range of surgical procedures. The Candy Cane Stirrups have adjustable heights and angles for gynecologic, urologic or laparoscopic procedures effectively. These SchureMed stirrups have a 350-pound weight capacity.

### Features

- Adjust vertically from 28" to 44" (71cm-112cm) with an 8" (20cm) leg clearance.
- Stirrups height can be locked by simply twisting the ergonomic handle.
- Connect surgical table stirrups to OR/surgical table side rail with Schure Socket XL (sold separately).

