



The CareFusion Alaris 8120 PCA Module works with the Alaris 8015 PC unit. The Alaris 8120 Module is used for delivery of pain management medications via patient-controlled dosing in a locked clear case.

The Alaris Medley 8120 PCA Syringe module is part of a comprehensive IV medication safety system with a common user interface shared across multiple infusion modalities. The 8120 is designed to work with the Alaris PC Unit as well as with other Alaris System modules.

The Alaris® PCA module is a syringe-based PCA infusion device that provides quick access to total dose history and monitors trending data giving you additional tools to analyze and address a patient's response to pain therapy.

Use this module to customize the Alaris 8000 or 8015 PC Unit Infusion Pumps

Features

- Reduced respiratory failure
- Decreased adverse drug events
- Monitoring for all at-risk patients



Specifications

Dimensions	Height: 15 in Width: 4.5 in Depth: 7.5 in Weight: 5.5lbs
Bolus Volume After Occlusion, Maximum	High: .997mL Low: .396 mL
Bolus Delivery Rate	150 mL/hr 75-500 mL/hr (limited by syringe size)
Lockout Interval	1 - 180 munte in 1 minute increments Min/max 1 - 180 minutes
Max. Accumulated Dose Range	1 hour limit Disabled/1, 2 or 4 hour limit
Max Rate	999 mL/hr 0.1 - 99.9 mL/h in 0.1 mL/h increments; 100 - 999 mL/h in 1 mL/h increments
Occlusion Pressure Set Point	High: 800 mmHg Low: 200 mmHg Medium: 500 mmhm

