



The 300 watt Xenon light source with a special four-channel turret that accommodates Olympus, Storz, ACMI/ Circon, and Wolf cable connectors. Additional features include intensity control and removable modular lamp. 300 watt xenon bulb has 188 hours on it. 120 VAC. ITEM is Used. Fully tested and functional



# Specifications

## General

Illumination level setting: Continuous

Mode: manual

Power needed, vac: 115/230

Hz:50/60

Dimensions: 14.5x 25.4 x 30.5 (5.7x 10 x 12)

Type: Multipurpose

Equivalent light cables: ACMI, Olympus, Storz, Wolf

High-intensity light: Yes

Lamp used: Xenon - for details please see comments.

Watts: 300

# of lamps: 1

Color temp, ÅK: 6000

Lamp duration, hr: 500

Other attributes: Can be used with any video camera; digital intensity read-out; universal turret. Meets CSA, IEC, and UL requirements