



The **ConMed AirSeal iFS** is an insufflation management system that can provide stable pneumoperitoneum, smoke evacuation, and valve-free access to the abdominal cavity during robotic and laparoscopy. This insufflator enables surgeons to operate at lower pressures or use standard high-flow insufflation. The AirSeal iFS lets the user select 3 different modes of operation that include: AirSeal Mode, Smoke Evacuation Mode, and Standard Insufflation Mode.

Features

- Stable Pneumoperitoneum, Smoke evacuation, and High-Flow Insufflation.
- Valve-free access to the abdominal cavity.
- Smoke evacuation and filtration with 0.01 μ ULPA Filter.
- Can use with conventional trocars.
- Can be used on adult and pediatric patients \geq 20kg.

