



The Beamer CE600 was specifically designed for endoscopic use and offers multiple therapies from a single platform. Combining sophisticated technology with unrivalled ease of use, the Beamer CE600 assists in the delivery of outstanding clinical results whatever your level of expertise in electrosurgery. Integrating built-in intelligence, innovation and an intuitive user interface, Beamer CE600 allows the quick selection of desired clinical effect, and delivers consistent repeatable clinical results.

Features

- Clear display; focuses on essentials
- Easy to operate; quick step control
- ABC® option
- Irrigation pump option
- Convenient storage
- Hands-free double foot pedal option



Specifications

Dimensions	(W x H x D) 390 x 182 x 435mm
Weight	8.3kg
Power requirements	100–240 V ₋ + 10%
Mains currents	max. 6.3 A
Mains fuse	6.3 A
Nominal frequencies (mains)	50 / 60 Hz
Power input standby mode	40 VA
Max. power input	600 VA
Output Data	
Cutting power(pulse)	max. 350 watts
Coagulation power	max. 120 watts
HF Frequencies	350 - 850 kHz
Classification acc. to MDD	IIb
Type of applied part	CF; defibrillation-proof
Approval/mark of conformity	CE 0297, in compliance with 93/42/EEC
LF leakage currents	complying to IEC 601, Part 1 (testing in connection with Beamer System)
Duty type	INT 10 s/30 s (= duty factor of 25%)