



CONMED's ClearView™ Smoke Evacuation System is a compact, quiet and effective smoke evacuation solution with state-of-the-art ULPA filtration to combat airborne hazards by removing 99.999% of the smoke contaminants. Surgical smoke contains over 80 chemicals that have the same level of mutagenicity as 27-30 cigarettes.

## Features

- Reduces noise pollution in the Operating Room.
- Compact, easy-to-integrated design allows for easy placement on shelf, counter, or stacking.
- System only activates when electrocautery is used for seamless integration with System 5000 ESU and other ESUs.
- Features three flow settings (Low/Lap, Medium, and High/Turbo) to meet the needs of different procedure types.
- Three port filter design accommodates a variety of tubing sizes, providing versatility in procedure and accessory setup.
- Up to 35 hour filter life.
- Delivers unsurpassed value by automatically tracking and recording filter life.
- Advanced, four stage filtration uses pre-filter, ULPA filter, carbon filter, and expanded foam to protect against particulates and microorganisms.



Make your OR smoke-free with CONMED's ClearView™ Smoke Evacuation System.

### Effective Evacuation

System features 4-stage ULPA filtration to remove 99.999% of the smoke contaminants from the surgical site<sup>1</sup>.

### Quiet Operation

Reduces noise pollution in the OR. Compact, Easy-to-Integrate Design Allows for easy placement on shelf, counter or stacking.

### Automatic Activation

System only activates when electro-cautery is used for seamless integration with System 5000 ESU and other ESUs.

### Multiple Settings for Multiple Procedures

Features three flow settings (Low/Lap, Medium & High/Turbo) to meet the needs of different procedure types.

### 3-Port Filter Design

Accommodates a variety of tubing sizes, providing versatility in procedure and accessory setup.

### Up to 35 Hour Filter Life

Delivers unsurpassed value by automatically tracking and recording filter life.



## Specifications

### Dimensions

Size: 16" x 11" x 6" (40.6 cm x 27.9 cm x 15.2 cm).

Shipping size: 21" x 16" x 15" (53.3 cm x 40.6 cm x 38.1 cm).

Shipping Weight: 19 lbs (8.6 kg).

### Power requirements

120 V; 60 Hz.