



The **Covidien McGRATH® MAC** can be used as a direct or indirect view laryngoscope. It utilizes all new single-use McGRATH® MAC blades made from robust optical polymer. It's supported by an internal reinforced alloy chassis and CameraStick™ which delivers 'steel-like' rigidity. Its super-slim blade profile brings improved patient access and limits dental contact. Its unique, minute-by-minute countdown battery technology can last for weeks.

Features

- **Easy to clean:** The screen and handle are fully immersible for high-level disinfection, and compatible with the latest vaporized hydrogen peroxide sterilization systems.
- **Slim Blade:** The slim blade profile provides improved patient access. Blades are packaged sterile and are disposable.
- **Portable:** The compact, cable-free device is ready when you need it.
- **Durable:** Drop-tested to 2 meters, the core of the McGRATH® MAC enhanced direct laryngoscope is reinforced with a steel alloy chassis.
- **Cost-effective:** The surprisingly affordable system and blades make the McGRATH® MAC enhanced direct laryngoscope as accessible as your conventional laryngoscope, with the flexibility of video laryngoscopy.



Specifications

Physical Specifications

Size: 180mm x 68mm x 110mm

Weight: 200g

Power: Proprietary 3.6V Lithium Battery Pack (c. 250 minutes)

Light source: High intensity LED

Display: 2.5" LCD color display

Camera: CMOS

Materials: Durable medical grade thermoplastics with reinforced structural alloy core.

The device and packaging are latex free

