



The **Cryomedics KRmed FO 101** is a colposcope that assists in a visual exam of your cervix that your physician will likely perform if you receive an abnormal result from a pap smear. Abnormal results might include the presence of HPV (human papillomavirus), cervical dysplasia (precancerous abnormal cells), cervical cancer, or squamous cells. A pap smear is not enough to diagnose the abnormal cells; it is just a screening to detect them. Therefore, more testing is required to check for or rule out HPV, precancerous lesions, or cancer. The colposcopy is usually the first follow-up exam your physician will perform.

Features

- Adjustable-tilt optical head with adjustable focus and 20X, 14 mm eyepieces
- Light intensity dial with off, low and high
- Uses 20 V 150 W halogen bulb
- 20.5" horizontal arm and 8.5" descending arm
- Triangular base with casters
- 7-foot power cord with grounded plug
- Colposcope Fiberlight
- Base: 22.5" x 22.5".
- Height: 56.5" (minimum)

