



The **Datex Ohmeda D-Tec Plus** desflurane anesthesia vaporizer is designed to add Suprane (desflurane) inhalant anesthetic vapor to the medical gases supplied to a patient. The vaporizer must be supplied with dry medical gas and it must be connected to a suitable electrical mains supply. The built-in filler assembly incorporates an agent specific fitting which helps ensure that only a Suprane (desflurane) specific bottle equipped with a "Saf-T-Fill" valve can be inserted to fill the vaporizer.

## Features

- Weight (Empty): 10 kg / 22 lbs
- Dimensions (HxWxD): 25cm x 11cm x 23.5cm
- Sump Capacity: 425 ml
- Flow Range: 0.2 to 10 liters/minute
- Concentration Range: 1 to 18%
- Flow Resistance: Less than 50 cm H<sub>2</sub>O at 5 liters/minute; Less than 150 cm H<sub>2</sub>O at 15 liters/minute.

