



The **GE Datex Ohmeda Tec 5** Sevoflurane anesthesia vaporizers are designed for out of circuit use in continuous flow techniques of inhalation anesthesia. The Tec 5 is temperature, flow, and pressure compensated so that its output remains relatively constant despite cooling due to evaporation. The Tec 5 is designed to be used on Selectatec Series Mounted Manifolds.

## Features

- Weight: 7 kg / 15.43 lbs
- Dimensions (HxWxD): 25cm x 11cm x 16.4cm
- Liquid Capacity (Dry wick): 300 Milliliters
- Liquid Retained by Wick: 75 Milliliters
- Calibration: Vaporizers are calibrated at 21° C using an oxygen carrier gas at a flow of 5 liters/minute and they are temperature, flow and pressure compensated within the specified operating range.

