



The Drager Air-Shields Isolette C450 QT infant incubator has been developed to address the unique requirements of infants. It provides thermally controlled environment required for any neonate. The Isolette C450 features both baby and air modes of operation. The important addition of the Quiet Touch (QT) helps reduce the adverse stimulation shown to have an effect on the well-being of the neonates.

Features

- Oval shaped port doors for increased access.
- Removable microprocessor controller with keypad lock and membrane touchpad.
- Iris ports for ventilation tubing support Typical manufacturer's specifications and images.
- Quiet Touch(QT) combines proven technology with innovative features to minimize excessive stimulation to the infant while, at the same time, reducing stress to parents and attending staff.



Specifications

Dimensions

Height: 55 in. (140 cm)

Width: 38 in (97 cm)

Depth 20 in (51 cm)

Weight: 168 lbs (76 kg)

Physical Characteristics

*Note: Wheels can add 3 each (6 max) to width or depth. Unit detaches from stand.

Adjustable mattress tilts $\pm 9^\circ \pm 1^\circ$

Wheels: 5" diameter, Antistatic type.

Mattress

Width: 25 in (63 cm)

Depth: 13.75 in (35 cm)

Physical characteristics with guard rail option

Width: 45.25 in. (115 cm)

Depth: 22 in. (56 cm)

Weight: of guard rail 13.5 lbs. (6 kg)