

Features

-
- A mobile medical unit, likely a mobile X-ray or ultrasound machine, featuring a blue and white design with a large screen and a control panel. It is mounted on a four-wheeled base with casters. The unit has a large, adjustable screen and a control panel with various buttons and a small display. The overall design is sleek and professional, typical of medical equipment.

Specifications

Dimensions

Height: 74.4" to 87" (189.6 to 221 cm)

Width: 29.5" (75 cm)

Length: 12.4" (31.5 cm)

Weight: 293.2 lbs (133 kg)

Radiant Warmer

Heat level 3: 10 mW/cm²

Heat level 10: 30 mW/cm²

Lamps

Work light: 30 W

Night light: 9W

Skin Temperature Measurements

Sensor: Drager sensors

Measuring Range: 59°F to 107.6°F (15°C to 42°C)

Accuracy: ± 0.1°C

Setting Range: 94°F to 99.5°F (35°C to 37.5°C)

Mattress Heater (Babytherm 8010 only)

Temperature Measurement: 41°F to 113°F (5°C to 45°C)

Measuring range/display range accuracy: ± 0.5°C between 20 and 42°C

Setting Range: 86°F to 101.3°F (30°C to 38.5°C)

Phototherapy (option)

Halogen 12V / 50W

Opening Angle: 24°

Radiation intensity at a distance of 80 cm: 9.5 W/m²

Useful Area: 15." x 7.8" (40 x 20 cm)