



The **Dräger Vapor 2000 Sevoflurane** is an unheated, calibrated anesthesia vaporizer for the enrichment of dry, fresh medical gas in anesthesia delivery systems with accurately-controlled concentrations of Sevoflurane vapor. The Vapor 2000 is connected between the fresh-gas flow-control unit and the fresh-gas outlet.

## Features

- Weight (*Filled*): 8.8 kg / 19.4 lbs
- Dimensions (*HxWxD*): 24.6cm x 10.8cm x 18.8cm
- Liquid Capacity (*Dry Wick*): ~360 mL
- Liquid Capacity (*Moist Wick*): ~360 mL
- Flow Range: 0.25 to 15 L/min
- Maximum Angle of Tilt: 30°

