



Features

- Wide-angle lens design for a larger field-of-view
- Unparalleled brightness
- Sharp images throughout the entire field-of-view
- Lightweight construction
- Sapphire lenses for improved scratch resistance compared to conventional glass
- Compatible with Smith & Nephew autoclavable and non-autoclavable Camera Heads
- Validated for both steam and chemical (EtO and Sterrad™) sterilization modalities



Specifications

Focus Ring Blue

Focal Length 35 mm

Outer Diameter 4.0 mm

Direction of View 30°

Field Of View 115°

Working Length 160mm

Light Post Direction Opposite

Cannula Connection J-Lock

Recommended Light Cable 5.0mm

