



The **ERBE VIO 100 C** is an electrosurgical unit that fulfills all basic requirements of monopolar cutting and monopolar and bipolar coagulation. The C models are designed for outpatient clinics and practices. The VIO 100 C is mostly used on interventions performed on an outpatient basis in the fields of dermatology, plastic surgery, ENT, gynecology, casualty surgery, and ophthalmology.

## Features

- Auto Cut - Minimal necrosis and reproducible cutting quality.
- Dry Cut - Cut with initial intense hemostasis.
- Soft Coag - Coagulation without carbonization and minimal sticking.
- Force Coag - Rapid and effective coagulation.
- Bipolar Soft Coag - Low voltage bipolar coagulation with minimal sticking and carbonization.



