



The **Ecolab ORS - 1058HS** hush slush machine delivers velvet-soft slush during procedures that require it. The ORS-1058HS uses proprietary Auto-Stir technology to automatically stir the solution while cooling. This eliminates concerns with conventional slush that may have sharp edges which can cause trauma to delicate tissues.

## Features

- Auto-stir device that automatically stirs and creates soft surgical slush.
- Disc-Drape that automatically snaps onto the Auto-Stir device for optimal sterile coverage.
- Caster wheels allow for easy transportation that lock to minimize inadvertent movement.
- Quiet operation offering less disruption to the operating room staff.
- Optional Ecolab ORS Fluid Warming System.



## Specifications

### Dimensions

**Height:** 36" (91.44 cm)

**Width:** 18" (45.72 cm)

**Depth:** 18" (45.72 cm)

### Basin Capacity

5.8 Liters (Rectangle)

### Fluid Level

1 to 3 Liters

