



The **Ecolab ORS-2057D** solution warmer is a fluid warming system. It provides immediate access to warm irrigation fluid at a visible and controlled temperature to increase patient safety and reduce the risk of hypothermia. The Ecolab ORS-2057D temperature settings can be adjusted and locked at a range from room temperature to 120°F (49°C).

Features

- Optional Portable Stand.
- Basin Capacity: 5.7 Liters.
- Dimensions: 19.13" H x 5.25" W x 4" D (48.59 cm H x 13.34 cm W x 10.2 cm D).
- Minimum Fluid Level: 0.75 liters.

