



The **Ecolab ORS-2066D** solution warmer is a fluid warming system. It provides immediate access to warm irrigation fluid at a visible and controlled temperature to increase patient safety and reduce the risk of hypothermia. The Ecolab ORS-2066D temperature settings can be adjusted and locked at a range from room temperature to 120°F (49°C).

Features

- Optional Portable Stand.
- Basin Capacity: 6.6 Liters.
- Dimensions: 11.5" H x 9.25" W x 4" D (29.21 cm H x 23.5 cm W x 10.2 cm D).
- Minimum Fluid Level: 0.75 liters.

